



Anti-Vaspin antibody ab82138

画像数 1

製品の概要

製品名	Anti-Vaspin antibody
製品の詳細	Goat polyclonal to Vaspin
由来種	Goat
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	<p>Synthetic peptide corresponding to Human Vaspin aa 211-224 (internal sequence) (Cysteine residue).</p> <p>Sequence:</p> <p>RARWKHEFDPNVTK</p> <p>Database link: NP_776249.1</p> <p> Run BLAST with  Run BLAST with</p>
ポジティブ・コントロール	WB: Human adipose tissue lysate.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	<p>pH: 7.30</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituents: Tris buffered saline, 0.5% BSA</p>
精製度	Immunogen affinity purified
特記事項 (精製)	ab82138 is purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

ポリ/モノ ポリクローナル
アイソタイプ IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab82138の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.3 - 1 µg/ml. Predicted molecular weight: 47 kDa. 1 hour primary incubation is recommended for this product.

ターゲット情報

機能	May modulates insulin action conceivably only in the presence of its yet undefined target proteases in white adipose tissues.
組織特異性	Expressed in visceral adipose tissues.
配列類似性	Belongs to the serpin family.
細胞内局在	Secreted.

画像



Anti-Vaspin antibody (ab82138) at 0.5 µg/ml + Human adipose lysate (in RIPA buffer) at 35 µg

Predicted band size: 47 kDa

Observed band size: 47 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors