


Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] ab14734

KO 評価済

★★★★★ [52 Abreviews](#) [512 References](#) [画像数 12](#)

製品の概要

| | |
|--------------|--|
| 製品名 | Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] |
| 製品の詳細 | Mouse monoclonal [20B12AF2] to VDAC1/Porin + VDAC3 |
| 由来種 | Mouse |
| アプリケーション | 適用あり: WB, ICC/IF, Flow Cyt |
| 種交差性 | 交差種: Mouse, Rat, Cow, Human 交差が予測される動物種: Sheep, Goat, Cat, Dog, Pig, Drosophila melanogaster, Fish, Quail, Common marmoset, Dogfish, Catshark  |
| 免疫原 | Recombinant full length protein. This information is proprietary to Abcam and/or its suppliers. |
| ポジティブ・コントロール | WB: Isolated mitochondria from human, cow, rat and mouse heart. HepG2 cell lysate. ICC/IF: HeLa cells. Human fibroblasts. Flow Cyt: HepG2 cells. |
| 特記事項 | <p>This antibody clone [20B12AF2] is manufactured by Abcam.</p> <p>If you require this antibody in a different buffer formulation or a different conjugate for your experiments, please contact orders@abcam.com or you can find further information here.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|-------|--|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Store at +4°C. Do Not Freeze. |
| バッファー | pH: 7.50 Preservative: 0.02% Sodium azide Constituents: 0.36% HEPES, 0.88% Sodium chloride |
| 精製度 | IgG fraction |

| | |
|----------|---|
| 特記事項(精製) | Near homogeneity as judged by SDS-PAGE. The antibody was produced in vitro using hybridomas grown in serum-free medium, and then purified by biochemical fractionation. |
| ポリ/モノ | モノクローナル |
| クローン名 | 20B12AF2 |
| アイソタイプ | IgG2b |
| 軽鎖の種類 | kappa |

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab14734の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|------------|---|
| WB | ★★★★★ (18) | Use a concentration of 1 µg/ml. Detects a band of approximately 39 kDa. |
| ICC/IF | ★★★★★ (4) | Use at an assay dependent concentration. |
| Flow Cyt | | Use 1µg for 10 ⁶ cells. ab170192 - Mouse monoclonal IgG2b, is suitable for use as an isotype control with this antibody. |

ターゲット情報

細胞内局在 VDAC1/Porin: Mitochondrion outer membrane. Cell membrane. VDAC3: Mitochondrial outer membrane.

画像



Western blot - Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734)

All lanes : Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734) at 1/1000 dilution

Lane 1 : Wild-type Hap1 cell lysate

Lane 2 : VDAC1 knockout Hap1 cell lysate

Lane 3 : Wild-type HEK-293T cell lysate


Lane 4 : VDAC1 knockout HEK-293T cell lysate

Lysates/proteins at 20 µg per lane.

Lanes 1 -4: Merged signal (red and green). Green - ab14734 observed at 31 kDa. Red - loading control, **ab181602** observed at 37 kDa. The lower band very close to VDAC1 is likely to be

VDAC3.

ab14734 was shown to react with VDAC1/Porin + VDAC3 in wild-type HEK-293T cells. Loss of signal was observed when knockout cell line [ab255444](#) (knockout cell lysate [ab263839](#)) was used. Wild-type and VDAC1 / Porin knockout samples were subjected to SDS-PAGE. ab14734 and Anti-GAPDH antibody [EPR16891] - Loading Control ([ab181602](#)) were incubated overnight at 4°C at 1 in 1000 dilution and 1 in 20000 dilution respectively. Blots were developed with Goat anti-Mouse IgG H&L (IRDye® 800CW) preadsorbed and Goat anti-Rabbit IgG H&L (IRDye® 680RD) preadsorbed secondary antibodies at 1 in 20000 dilution for 1 hour at room temperature before imaging.




Immunocytochemistry/ Immunofluorescence - Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734)

This image is courtesy of an Abreview by Michiel Krols.

ab14734 staining VDAC1/Porin + VDAC3 in human HeLa (Human epithelial cell line from cervix adenocarcinoma) cells by ICC/IF (Immunocytochemistry/immunofluorescence).

Cells were fixed with paraformaldehyde, permeabilized with 0.1% Triton X-100 for 3 minutes and blocked with 0.2% serum for 60 minutes at 22°C. Samples were incubated with primary antibody (1/200 in 0.5% BSA and 0.02% Triton X100 in PBS) for 16 hours at 4°C. An FITC-conjugated goat anti-mouse IgG polyclonal was used as the secondary antibody at a dilution of 1/200.



Western blot - Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734)

All lanes : Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734) at 1/1000 dilution

Lane 1 : N-GST tagged Recombinant Human VDAC1 (aa1 to 283) protein with NFDm/TBST

Lane 2 : N-GST tagged Recombinant Human VDAC2 (aa 1 to 294) protein with NFDm/TBST

Lane 3 : N-GST tagged Recombinant Human VDAC3 (aa 1 to 283) protein with NFDm/TBST

Lysates/proteins at 0.01 µg per lane.

Blocking peptides at 5 % per lane.


Observed band size: 33 kDa

Exposure time: 10 seconds

Diluting buffer: 5% NFDM /TBST

N-GST tagged Recombinant Human VDAC1 protein is available as **ab132481**

N-GST tagged Recombinant Human VDAC2 protein is available as **ab152793**



Western blot - Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734)

All lanes : Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734) at 1/1000 dilution

Lane 1 : Wild type HAP1 whole cell lysate


Lane 2 : VDAC1 knockout HAP1 whole cell lysate

Lysates/proteins at 20 µg per lane.

Observed band size: 33 kDa

Exposure time: 5 seconds

Blocking and dilution buffer: 5% NFDM /TBST.



Western blot - Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734)

All lanes : Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734) at 1/1000 dilution

Lane 1 : HEK-293 (Human embryonic kidney epithelial cell) whole cell lysate with NFDM/TBST

Lane 2 : HAP1 (Wildtype control Human chronic myelogenous leukemia near-haploid cell line) whole cell lysate with NFDM/TBST

Lane 3 : MEF (Mus musculus Embryo Fibroblast) whole cell lysate with NFDM/TBST

Lane 4 : Rat brain tissue lysate with NFDM/TBST


Lysates/proteins at 20 µg per lane.

Blocking peptides at 5 % per lane.

Observed band size: 33 kDa

Exposure time: 5 seconds

Diluting buffer: 5% NFDM /TBST



Western blot - Anti-VDAC1/Porin + VDAC3 antibody
[20B12AF2] (ab14734)

Sun Y et al Voltage-dependent anion channels (VDACs) recruit Parkin to defective mitochondria to promote mitochondrial autophagy. J Biol Chem. 2012 Nov 23;287(48):40652-60. doi: 10.1074/jbc.M112.419721.


(A) The Anti-VDAC1/Porin + VDAC3 antibody recognizes VDAC1 and VDAC3 but not VDAC2.

Lysates from wild-type MEFs and VDAC1/3^{-/-} MEFs were analyzed by western blotting with Anti-VDAC1/Porin + VDAC3 (top panel) or anti Tom20 antibodies (bottom panel).

(B) The Anti-VDAC1/Porin + VDAC3 antibody does not recognize VDAC2 in VDAC1/3^{-/-} MEFs.

HA-tagged VDAC1, 2 or 3 were expressed in VDAC1/3^{-/-} MEFs, immunoprecipitated with anti-HA antibodies and analyzed by western blotting with anti-HA antibodies (top panel) or Anti-VDAC1/Porin + VDAC3 antibodies (bottom panel). The Anti-VDAC1/Porin + VDAC3 antibodies recognize HA-VDAC1 and HA-VDAC3, but not HA-VDAC2. Note that HA-VDAC3 migrates at a lower molecular weight than HA-VDAC1 and HA-VDAC2.

Molecular weights are indicated in kDa.



Western blot - Anti-VDAC1/Porin + VDAC3 antibody
[20B12AF2] (ab14734)


Li et al PLoS One. 2016 Feb 22;11(2):e0149728. doi: 10.1371/journal.pone.0149728. eCollection 2016. Fig 5.

Immunoblot analysis of endogenous VDAC1/Porin and VDAC3 in subcellular fractions of MA-10 cell lysates. VDAC1/Porin and VDAC3 are used as a mitochondrial marker proteins.

Lane 1: Whole cell lysate

Lane 2: Cytosolic fraction

Lane 3: Mitochondrial-enriched fraction



Western blot - Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734)

All lanes : Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734)

Lane 1 : Isolated mitochondria from human heart at 15 μ g

Lane 2 : Isolated mitochondria from bovine heart at 6 μ g


Lane 3 : Isolated mitochondria from rat heart at 30 μ g

Lane 4 : Isolated mitochondria from mouse heart at 30 μ g

Lane 5 : HepG2 (Human liver hepatocellular carcinoma cell line) cell lysate at 30 μ g

Observed band size: 37 kDa

Extra bands in the mouse sample (lane 4) are due to the reaction of the secondary antibody with residual mouse blood in the heart tissue, as it is very difficult to entirely remove the blood from these small organs.




Flow Cytometry - Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734)

Overlay histogram showing HepG2 (Human liver hepatocellular carcinoma cell line) cells stained with ab14734 (red line).

The cells were fixed with 4% paraformaldehyde (10 min) and then permeabilized with 0.1% PBS-Tween for 20 min. The cells were then incubated in 1x PBS / 10% normal goat serum / 0.3M glycine to block non-specific protein-protein interactions followed by the antibody (ab14734, 1 μ g/1x10⁶ cells) for 30 min at 22°C. The secondary antibody used was DyLight® 488 goat anti-mouse IgG (H+L) ([ab96879](#)) at 1/500 dilution for 30 min at 22°C. Isotype control antibody (black line) was mouse IgG2b [PLPV219] ([ab91366](#), 2 μ g/1x10⁶ cells) used under the same conditions.


Acquisition of >5,000 events was performed.

This antibody gave a positive signal in HepG2 cells fixed with 80% methanol (5 min)/permeabilized in 0.1% PBS-Tween used under the same conditions.



Immunocytochemistry/ Immunofluorescence - Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734)

Immunofluorescence using ab14734 at 0.2 µg/ml on human fibroblasts (red).
Nuclei were labeled with DAPI (blue).



Western blot - Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734)
Image courtesy of an anonymous abreview.

Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734) at 1/500 dilution + MCF7 (Human breast adenocarcinoma cell line) cell lysate at 10 µg

Secondary


HRP-conjugated goat anti-mouse polyclonal at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Exposure time: 3 minutes

0.5% TBS-tween + Lait 5% NaN₃ for 16 hours at 4°C.



Western blot - Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734)
This image is courtesy of an anonymous Abreview

All lanes : Anti-VDAC1/Porin + VDAC3 antibody [20B12AF2] (ab14734) at 1/5000 dilution

Lanes 1-2 : Rat brain cell lysate (homogenate)

Lanes 3-4 : Rat brain cell lysate (mitochondrial)

Lysates/proteins at 30 µg per lane.

Secondary

All lanes : HRP conjugated sheep anti-mouse IgG

Developed using the ECL technique.

Performed under reducing conditions.

Observed band size: 39 kDa

Exposure time: 5 minutes

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors