

### Anti-VENTX antibody ab124523

面格新 1

医薬用外劇物

#### 製品の概要

製品名	Anti-VENTX antibody
製品の詳細	Rabbit polyclonal to VENTX
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Mouse, Rat, Human
免疫原	Synthetic peptide surrounding amino acid 76 of Human VENTX (UniProt O95231).
ポジティブ・コントロール	Jurkat cell lysate, 3T3 mouse lysate and Rat kidney tissue lysate.
特記事項	<p>This product is manufactured by BioVision, an Abcam company and was previously called 3286 VENTX Antibody. 3286-100 is the same size as the 100 µg size of ab124523.</p> <p>The antibody solution should be gently mixed before use.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: 69% PBS, 30% Glycerol (glycerin, glycerine), 0.5% BSA, 0.15% EDTA
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション

### The Abpromise guarantee

**Abpromise保証は、**次のテスト済みアプリケーションにおけるab124523の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.5 - 4 µg/ml. Predicted molecular weight: 27 kDa.

## ターゲット情報

機能	May be involved in ventralization.
組織特異性	Expressed in bone marrow of patients recovering from chemotherapy. Also expressed in an erythroleukemia cell line.
配列類似性	Contains 1 homeobox DNA-binding domain.
細胞内局在	Nucleus.

## 画像



**All lanes :** Anti-VENTX antibody (ab124523) at 4 µg/ml

**Lanes 1-2 :** Jurkat cell lysate

**Lane 3 :** 3T3 cell lysate

**Lane 4 :** Rat kidney tissue lysate

**Predicted band size:** 27 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors