


Anti-WDR36 antibody ab73548

★★★★☆ [1 Abreviews](#) [1 References](#)

製品の概要

| | |
|----------|---|
| 製品名 | Anti-WDR36 antibody |
| 製品の詳細 | Rabbit polyclonal to WDR36 |
| 由来種 | Rabbit |
| アプリケーション | 適用あり: IHC-Fr, IHC-P, ICC/IF, WB |
| 種交差性 | 交差種: Human 交差が予測される動物種: Mouse  |
| 免疫原 | Synthetic peptide derived from N terminal domain of human WDR36 |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | Constituent: Whole serum |
| 精製度 | Whole antiserum |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab73548の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|--|
| IHC-Fr | | Use at an assay dependent concentration. |
| IHC-P | | Use at an assay dependent concentration. |
| ICC/IF | ★★★★☆ (1) | Use at an assay dependent concentration. |
| WB | | Use at an assay dependent concentration. |

ターゲット情報

| | |
|--------------|--|
| 関連性 | WDR36 is a member of the WD repeat protein family. WD repeats are minimally conserved regions of approximately 40 amino acids typically bracketed by gly-his and trp-asp (GH-WD), which may facilitate formation of heterotrimeric or multiprotein complexes. WDR36 is involved in T cell activation and is highly co-regulated with IL2. Mutations in WDR36 are associated with adult-onset primary open-angle glaucoma (POAG). |
| 細胞内局在 | Nuclear |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors