

Anti-White Spot Syndrome Virus VP664 antibody ab26811

画像数 1

製品の概要

製品名	Anti-White Spot Syndrome Virus VP664 antibody
製品の詳細	Rabbit polyclonal to White Spot Syndrome Virus VP664
由来種	Rabbit
アプリケーション	適用あり: ELISA
種交差性	交差種: White spot syndrome virus
免疫原	Synthetic peptide located on the putative domain of White spot syndrome virus VP664.
ポジティブ・コントロール	synthetic WSSV VP664 peptideantigen.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab26811の使用に適用されず**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ELISA		1/1000.

ターゲット情報

関連性 White spot syndrome virus [WSSV] is one of the most devastating shrimp pathogens, and it has caused serious damage to the shrimp culture industry. This virus can infect several crustacean species, including shrimp, crab, and crayfish. VP664, is one of the viral structural proteins. The existence of the corresponding 18-kb transcript was confirmed by sequencing analysis of reverse transcription-PCR products, which also showed that vp664 was intron-less.

画像



ELISA using ab26811 at varying antibody concentrations and antigen concentration at 1000 ng/mL.

ELISA - Anti-White Spot Syndrome Virus VP664
antibody (ab26811)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors