

Anti-Wnt1 antibody ab63934

6 References [画像数 1](#)

医薬用外劇物

製品の概要

製品名	Anti-Wnt1 antibody
製品の詳細	Rabbit polyclonal to Wnt1
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide surrounding amino acid 31 of human Wnt1
ポジティブ・コントロール	WB: Jurkat cell lysate.
特記事項	<p>This product is manufactured by BioVision, an Abcam company and was previously called 3571 Anti-Wnt-1 Antibody. 3571-100 is the same size as the 100 µg size of ab63934.</p> <p>Wnt1 antibody has been purified by Protein A/G combination.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.01% Thimerosal (merthiolate) Constituents: PBS, 30% Glycerol (glycerin, glycerine), 0.5% BSA
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab63934の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 0.5 - 4 µg/ml. Detects a band of approximately 41 kDa (predicted molecular weight: 41 kDa).

ターゲット情報

機能	Ligand for members of the frizzled family of seven transmembrane receptors. In some developmental processes, is also a ligand for the coreceptor RYK, thus triggering Wnt signaling. Probable developmental protein. May be a signaling molecule important in CNS development. Is likely to signal over only few cell diameters.
配列類似性	Belongs to the Wnt family.
細胞内局在	Secreted > extracellular space > extracellular matrix.

画像



Western blot - Anti-Wnt1 antibody (ab63934)

All lanes : Anti-Wnt1 antibody (ab63934) at 1 µg/ml

All lanes : Jurkat cell lysate

Lysates/proteins at 35 µg per lane.

Predicted band size: 41 kDa

Observed band size: 41 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

-
- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors