


### Anti-ZBTB1 antibody ab115170

画像数 1

#### 製品の概要

---

製品名	Anti-ZBTB1 antibody
製品の詳細	Rabbit polyclonal to ZBTB1
由来種	Rabbit
アプリケーション	<b>適用あり:</b> WB
種交差性	<b>交差種:</b> Mouse, Human <b>交差が予測される動物種:</b> Rat, Rabbit, Chicken, Guinea pig, Cow, Dog, Turkey, Pig, Chimpanzee, Rhesus monkey, Gorilla 
免疫原	Synthetic peptide, corresponding to a region within amino acids 663 - 713 of Human ZBTB1 (Q942K1).
ポジティブ・コントロール	HeLa, Jurkat, 293T and NIH3T3 cell lysate
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

#### 製品の特性

---

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7 Preservative: 0.09% Sodium azide Constituent: 99% Tris citrate/phosphate
精製度	pH 7 to 8. Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

## アプリケーション


**The Abpromise guarantee**      **Abpromise保証は、** 次のテスト済みアプリケーションにおけるab115170の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/1000 - 1/5000. Predicted molecular weight: 82 kDa.

## ターゲット情報

機能	May be involved in transcriptional regulation.
配列類似性	Contains 1 BTB (POZ) domain. Contains 8 C2H2-type zinc fingers.
細胞内局在	Nucleus.

## 画像



**All lanes :** Anti-ZBTB1 antibody (ab115170) at 0.4 µg/ml

**Lane 1 :** HeLa cell lysate at 50 µg/ml

**Lane 2 :** HeLa cell lysate at 15 µg/ml

**Lane 3 :** Jurkat cell lysate at 50 µg/ml

**Lane 4 :** 293T cell lysate at 50 µg/ml

**Lane 5 :** Mouse NIH3T3 cell lysate at 50 µg/ml

**Predicted band size:** 82 kDa

**Exposure time:** 3 minutes

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors