

Anti-Zebrafish Basement Membrane From Intestinal Epithelial antibody [5F11] ab71312

1 References

製品の概要

製品名	Anti-Zebrafish Basement Membrane From Intestinal Epithelial antibody [5F11]
製品の詳細	Mouse monoclonal [5F11] to Zebrafish Basement Membrane From Intestinal Epithelial
由来種	Mouse
特異性	The transparency of the juvenile zebrafish and its genetic advantages make it an attractive model for the study of intestinal differentiation and renewal. 5F11 is a marker for the basement membrane of the intestinal epithelium.
アプリケーション	適用あり: ICC/IF, WB, IP
種交差性	交差種: Zebrafish
免疫原	Tissue, cells or virus corresponding to Zebrafish Zebrafish Basement Membrane From Intestinal Epithelial. Zebrafish Intestine.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: 0.02% Sodium azide Constituent: 99.98% PBS
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	5F11
アイソタイプ	IgG1

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab71312の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		
WB		
IP		

追加情報

ICC/IF: Use at an assay dependent dilution.
IP: Use at an assay dependent dilution.
WB: Use at an assay dependent dilution.
This antibody is also suitable for use in histology.

Not yet tested in other applications.
Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors