

Anti-ZGRF1 antibody ab122126

画像数 1

製品の概要

製品名	Anti-ZGRF1 antibody
製品の詳細	Rabbit polyclonal to ZGRF1
由来種	Rabbit
アプリケーション	適用あり: IHC-P
種交差性	交差種: Human
免疫原	Recombinant fragment corresponding to Human ZGRF1 aa 1741-1828. Sequence: ESEQLKELHALMKEDLTPTERVYVRKSIEQHKLGTRNRTLLKQ VRVVGVTCT AACPFPCMNDLKFPVVVLDECSQITEPASLLPIARFEC Run BLAST with Run BLAST with
ポジティブ・コントロール	Human placenta tissue
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.20 Preservative: 0.02% Sodium azide Constituents: 40% Glycerol (glycerin, glycerine), 59% PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors