

### CSNK2A1 peptide ab204855

#### 製品の詳細

製品名	CSNK2A1 peptide
精製度	> 97 % HPLC.
Animal free	No
由来	Synthetic
配列	RRRADDSDDDDD

#### 特性

Our **Abpromise guarantee** covers the use of **ab204855** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

製品の状態	Lyophilized
備考	ab204855 (CKII peptide) can be utilized as a substrate for <b>ab105125</b> (Active human CSNK2A1 prime polypeptide full length protein) This product was previously labelled as CKII alpha

#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Avoid freeze / thaw cycle.
再構成	Dilute peptide in 20mM Tris-HCl, pH 7.5 solution to a final concentration of 1mg/ml

#### 関連情報

機能	Casein kinases are operationally defined by their preferential utilization of acidic proteins such as caseins as substrates. The alpha and alpha' chains contain the catalytic site. Participates in Wnt signaling. CK2 phosphorylates 'Ser-392' of p53/TP53 following UV irradiation.
配列類似性	Belongs to the protein kinase superfamily. Ser/Thr protein kinase family. CK2 subfamily. Contains 1 protein kinase domain.

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors