

### Human Atrial natriuretic factor peptide ab42331

#### 製品の詳細

製品名	Human Atrial natriuretic factor peptide
精製度	> 95 % HPLC.
アクセッション番号	<b><u>P01160</u></b>
Animal free	No
由来	Synthetic
生物種	Human
配列	SLRRSSCFGGRMDRIGAQSGLGCNSFRY
領域	124 to 151

#### 特性

Our **Abpromise guarantee** covers the use of **ab42331** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Blocking
製品の状態	Lyophilized

#### 前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Constituent: 0.1% Acetic acid
再構成	Resuspend in 0.1% acetic acid for required concentration

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors