

Human proBNP peptide ab42238

製品の詳細

製品名	Human proBNP peptide
精製度	> 95 % HPLC.
Animal free	No
由来	Synthetic
生物種	Human
配列	QVEQTSLEPLQESPRPTGVWKSREVA
領域	32 to 57

特性

Our **Abpromise guarantee** covers the use of **ab42238** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Blocking
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Constituent: 0.1% Acetic acid
再構成	Resuspend in 0.1% acetic acid for required concentration (vial content is 0.1 mg).

関連情報

機能	Cardiac hormone which may function as a paracrine antifibrotic factor in the heart. Also plays a key role in cardiovascular homeostasis through natriuresis, diuresis, vasorelaxation, and inhibition of renin and aldosterone secretion. Specifically binds and stimulates the cGMP production of the NPR1 receptor. Binds the clearance receptor NPR3.
組織特異性	Brain and also in atria, but at much lower levels than ANP.
配列類似性	Belongs to the natriuretic peptide family.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors