

Human Thymosin alpha 1 peptide ab42247

2 References

製品の詳細

製品名	Human Thymosin alpha 1 peptide
精製度	> 95 % HPLC.
Animal free	No
由来	Synthetic
生物種	Human
配列	SDAAVDTSSSEITTKDLKEKKEVVEEAEN
領域	2 to 29

特性

Our **Abpromise guarantee** covers the use of **ab42247** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Blocking - Blocking peptide for Anti-Thymosin alpha 1 antibody (ab14332)
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles. Constituent: 0.1% Acetic acid
再構成	Store at 2-8°C lyophilised, resuspend in 0.1% acetic acid.

関連情報

機能	Prothymosin alpha may mediate immune function by conferring resistance to certain opportunistic infections.
配列類似性	Belongs to the pro/parathymosin family.
翻訳後修飾	Covalently linked to a small RNA of about 20 nucleotides.
細胞内局在	Nucleus.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors