

Native Chicken IgY protein ab119138

1 References

製品の詳細

製品名	Native Chicken IgY protein
精製度	>= 97 % n/a. IgY fraction. Immunoglobulin fraction purified from egg yolk from eggs laid by hens that have not been immunized. Product is 0.2 µm filtered.
発現系	Native
タンパク質長	Full length protein
Animal free	No
由来	Native
生物種	Chicken

特性

Our [Abpromise guarantee](#) covers the use of **ab119138** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

製品の状態	Liquid
-------	--------

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C. Preservative: 0.1% Sodium azide Constituents: 0.14% Sodium phosphate, 0.88% Sodium chloride
------------	--

関連情報

関連性	Secreted immunoglobulins have two separable functions: one is to bind specifically to pathogenic molecules, the other is to recruit cells and molecules to destroy the opsonized pathogen. IgY has a tetrameric structure consisting of two identical light chains (Mr = 24 kDa) and two identical heavy chains (Mr = 73 kDa).
細胞内局在	Secreted

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors