

Native Clostridium difficile Toxoid B protein ab124002

画像数 1

製品の詳細

製品名	Native Clostridium difficile Toxoid B protein
精製度	> 95 % SDS-PAGE.
発現系	Native
タンパク質長	Full length protein
Animal free	No
由来	Native
予測される分子量	270 kDa

特性

Our **Abpromise guarantee** covers the use of **ab124002** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE ELISA
----------	-------------------

製品の状態	Lyophilized
-------	-------------

備考	Inactivation: Inactivated using formaldehyde. No cell rounding in the sensitive verocytotoxicity assay at 2µg/ml compared to active toxin lethal dose of 4pg/ml. Antigenicity: Maintained to be no less than 90% of original toxin strength as measured by direct ELISA.
----	---

This protein is derived from a pathogenic organism, and may be involved in a disease process, consequently exposure may have adverse health effects.

前処理および保存


保存方法および安定性	Shipped at 4°C. Store at +4°C. Avoid freeze / thaw cycle. Reconstitute for long term storage. pH: 7.00 Constituents: 1.19% HEPES, 5% Sucrose, 0.88% Sodium chloride
------------	---

再構成	Reconstitute in 250µl sterile distilled water.
-----	--

関連情報

関連性 Clostridium difficile, or C. difficile, a gram positive anaerobic bacterium, is the major causative agent of colitis and diarrhoea that can occur following antibiotic intake. C. difficile infection is one of the most common hospital acquired infections. When C. difficile colonize the gut, they release two potent toxins, toxin A and toxin B, which bind to certain receptors in the lining of the colon and ultimately cause diarrhoea and inflammation of the large intestine.

画像



SDS-PAGE showing ab124002 at 270kDa.

SDS-PAGE - Native Clostridium difficile Toxoid B protein (ab124002)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors