

Native Cow Bovine Serum Albumin protein ab64009

5 References

製品の詳細

製品名	Native Cow Bovine Serum Albumin protein
生理活性	This product contains:- Moisture: 1.4% Ash: 1.4% Heavy Metals: <0.44ppm Bovine IgG: Not detected Protein: 100% (Biuret) Albumin: 100% (electrophoresis)
精製度	> 96 % SDS-PAGE. Proprietary heat-shock process.
発現系	Native
アクセッション番号	<u>P02769</u>
タンパク質長	Full length protein
Animal free	No
由来	Native
生物種	Cow
配列	DTHKSEIAHRFKDLGEEHFKGLVLIAFSQYLQQCPFDEHVKL VNELTEFA KTCVADESHAGCEKSLHTLFGDELCKVASLRETYGDMADCCCE KQEPERNE CFLSHKDDSPDLPKLKPDPNTLCDEFKADEKKFWGKYLYEIA RRHPYFYA PELLYYANKYNGVFQECQAEDKGACLLPKIETMREKVLASS ARQLRCA SIQKFGERALKAWSVARLSQKFPKAEFVEVTKLVTDLTKVHK ECCHGDLL ECADDRADLAKYICDNQDTISSKLKECCDKPLLEKSHCIAEV EKDAIPEN LPPLTADFAEDKDVCKNYQEAKDAFLGSFLYEYSRRHPEYAV SVLLRLAK EYEATLEECCAADDPHACYSTVFDKCLKHLVDEPQNLIKQNC QFEKLGEY GFQNALIVRYTRKVPQVSTPTLVEVSRSLGKVGTRCCTKPES ERMPCTED YLSLILNRLCVLHEKTPVSEKVTCCTESLVNRRPCFSALTP DETYVPKA FDEKLFTFHADICTLPDTEKQIKKQTALVELLKHKPKATEEQ LKTVMENF VAFVDKCCAADDKEACFAVEGPKLVVSTQTALA
領域	25 to 607

配列の追加情報

Source = bovine serum

特性

Our **Abpromise guarantee** covers the use of **ab64009** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション

Western blot

Immunohistochemistry

製品の状態

Lyophilized

備考

The product is suitable to use in blocking experiments. This product contains:- Moisture: 1.4%
Ash: 1.4% Heavy Metals: <0.44ppm Bovine IgG: Not detected Protein: 100% (Biuret) Albumin:
100% (electrophoresis)

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

再構成

Reconstitute using deionized water.

関連情報

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors