

Native Herpes simplex virus HSV2 protein ab124006

製品の詳細

製品名	Native Herpes simplex virus HSV2 protein
発現系	Native
タンパク質長	Full length protein
Animal free	No
由来	Native
生物種	Herpes simplex virus

特性

Our **Abpromise guarantee** covers the use of **ab124006** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

製品の状態	Liquid
備考	<p>HSV2 is propagated in mammalian cells using a protocol optimised for maximum particle yield. At the appropriate time, infected cells are washed in saline then disrupted to release virus particles by freeze thaw and sonication.</p> <p>Inactivation: 0.1% Triton x-100.</p> <p>This protein is derived from a pathogenic organism, and may be involved in a disease process, consequently exposure may have adverse health effects.</p> <p>Abcam has not and does not intend to apply for the REACH Authorisation of customers' uses of products that contain European Authorisation list (Annex XIV) substances.</p> <p>It is the responsibility of our customers to check the necessity of application of REACH Authorisation, and any other relevant authorisations, for their intended uses.</p>

前処理および保存

保存方法および安定性	<p>Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.</p> <p>pH: 7</p> <p>Constituents: 99% PBS, 0.1% Triton-X-100</p>
------------	---

関連情報

関連性

Herpes simplex type 2 (HSV2) belongs to a family that includes HSV1, Epstein-Barr virus (EBV) and Varicella zoster (chicken pox) virus. HSV1 and HSV2 are extremely difficult to distinguish from each other. These viruses have a DNA genome, an icosahedral protein coat and are encased in a lipid membrane derived from the nuclear membrane of the last host. These viruses are capable of entering a latent phase where the host shows no visible sign of infection and levels of infectious agent become very low. During the latent phase the viral DNA is integrated into the genome of the host cell.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors