

Native Human IgG Fab fragment protein ab90352

2 References

製品の詳細

| | |
|-------------|---------------------------------------|
| 製品名 | Native Human IgG Fab fragment protein |
| 精製度 | > 95 % SDS-PAGE. |
| 発現系 | Native |
| タンパク質長 | Full length protein |
| Animal free | No |
| 由来 | Native |
| 予測される分子量 | 50 kDa |

特性

Our **Abpromise guarantee** covers the use of **ab90352** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|--|
| アプリケーション | SDS-PAGE |
| 製品の状態 | Liquid |
| 備考 | Protein Determination: Extinction Coefficient (E) 0.1% at 280nm = 1.35. Prepared from plasma shown to be non reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA approved tests. |

前処理および保存

| | |
|------------|---|
| 保存方法および安定性 | Shipped on Dry Ice. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 0.0536% PBS, 0.87% Sodium chloride |
|------------|---|

関連情報

| | |
|-----|--|
| 関連性 | Human IgG consists of four subclasses (1-4) that can be recognised by antigen differences in their heavy chains. Each subclass has different biological and physiochemical properties. The IgG |
|-----|--|

subclass may be preferentially produced in response to different antigens. The Fab fragment (fragment antigen binding) is the antigen region of the IgG molecule. It is composed of one constant and one variable domain of each of the heavy and the light chain.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors