

# Native Human IgG protein ab91102

8 References 画像数 1

医薬用外毒物

### 製品の詳細

製品名	Native Human IgG protein
精製度	> 95 % SDS-PAGE.
発現系	Native
タンパク質長	Full length protein
Animal free	No
由来	Native
生物種	Human
予測される分子量	150 kDa

### 特性

Our **Abpromise guarantee** covers the use of **ab91102** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Lyophilized
備考	Protein Determination: Extinction Coefficient (E) 0.1% at 280nm, 1cm pathway = 1.36.  Prepared from plasma shown to be non reactive for HBsAg, anti-HCV, anti-HBc, and negative for anti-HIV 1 & 2 by FDA approved tests.

### 前処理および保存


保存方法および安定性	Shipped at 4°C. Store at -20°C.  pH: 7.40 Preservative: 0.05% Sodium azide Constituents: 0.87% Sodium chloride, Phosphate Buffer
再構成	Centrifuge briefly the vial prior to opening. Reconstitute with sterile deionized water to a concentration of 47.6 mg/mL, followed by gentle swirling and/or vortexing to ensure adequate homogenization. If further dilution is required, use 20 mM Phosphate buffer, pH 7.4, with 150 mM NaCl and 0.05% sodium azide.

## 関連情報

細胞内局在

Secreted

## 画像



SDS-PAGE: 4-12% Bis-Tris NuPAGE gel

Lane 1. Molecular weight markers

Lane 2. 5 µg IgG (reduced/heated)

Lane 3. 10 µg IgG (reduced/heated)

Lane 4. 20 µg IgG (reduced/heated)

Lane 5. 5 µg IgG (non-reduced/no heat)

Lane 6. 10 µg IgG non-reduced/no heat)

Lane 7. 20 µg IgG non-reduced/no heat)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors