

Native Human Lactoferrin protein ab90354

製品の詳細

| | |
|-------------|----------------------------------|
| 製品名 | Native Human Lactoferrin protein |
| 精製度 | > 95 % SDS-PAGE. |
| 発現系 | Native |
| タンパク質長 | Full length protein |
| Animal free | No |
| 由来 | Native |
| 生物種 | Human |

特性

Our **Abpromise guarantee** covers the use of **ab90354** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|----------|----------|
| アプリケーション | SDS-PAGE |
| 製品の状態 | Liquid |

前処理および保存

| | |
|------------|---|
| 保存方法および安定性 | Shipped on Dry Ice. Store at -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 1.58% Tris HCl, 1.74% Sodium chloride |
|------------|---|

関連情報

| | |
|----|---|
| 機能 | <p>Transferrins are iron binding transport proteins which can bind two Fe(3+) ions in association with the binding of an anion, usually bicarbonate.</p> <p>Lactotransferrin has antimicrobial activity which depends on the extracellular cation concentration. Lactoferroxins A, B and C have opioid antagonist activity. Lactoferroxin A shows preference for mu-receptors, while lactoferroxin B and C have somewhat higher degrees of preference for kappa-receptors than for mu-receptors.</p> <p>The lactotransferrin transferrin-like domain 1 functions as a serine protease of the peptidase S60 family that cuts arginine rich regions. This function contributes to the antimicrobial activity.</p> |
|----|---|

| | |
|-------|--|
| 配列類似性 | Belongs to the transferrin family. Contains 2 transferrin-like domains. |
| 細胞内局在 | Secreted. |

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors