

Native pig Elastase protein ab95111

製品の詳細

製品名	Native pig Elastase protein
生理活性	ab95111 cleaves PAI-1 between the P3 and P4 residues to create a loop-inserted species that is non-inhibitory towards the target enzymes uPA and tPA. It can be used to cleave Human PAI-1 in any application requiring elastase cleavage.
精製度	> 95 % SDS-PAGE.
発現系	Native
タンパク質長	Full length protein
Animal free	No
由来	Native
生物種	Pig

特性

Our **Abpromise guarantee** covers the use of **ab95111** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Functional Studies
製品の状態	Liquid
備考	ab95111 cleaves PAI-1 between the P3 and P4 residues to create a loop-inserted species that is non-inhibitory towards the target enzymes uPA and tPA. It can be used to cleave Human PAI-1 in any application requiring elastase cleavage.

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C. pH: 7.40 Preservative: 0.02% Sodium azide Constituents: 2.38% HEPES, 0.58% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--

関連情報

機能	Efficient protease with alanine specificity but only little elastolytic activity.
配列類似性	Belongs to the peptidase S1 family. Elastase subfamily. Contains 1 peptidase S1 domain.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors