

Recombinant *E. coli* Beta-lactamase protein ab104926

画像数 1

製品の詳細

製品名	Recombinant <i>E. coli</i> Beta-lactamase protein
精製度	> 95 % SDS-PAGE. ab104926 was purified using conventional chromatography.
発現系	<i>Escherichia coli</i>
アクセッション番号	<u>P00811</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	<i>Escherichia coli</i>
配列	<pre>MGSSHHHHH SSGLVPRGSH MAPQQINDIV HRTITPLIEQ QKIPGMAVAV IYQGKPYFT WGYADIAKKQ PVTQQTLFEL GSVSKTFTGV LGGDAIARGE IKLSDPTTKY WPELTAKQWN GITLLHLATY TAGGLPLQVP DEVKSSDLL RFYQNWQPAW APGTQRLYAN SSIGLFGALA VKPSGLSFEQ AMQTRVFQPL KLNHTWINVP PAEEKNYAWG YREGKAVHVS PGALDAEAYG VKSTIEDMAR WVQSNLKPLD INEKTLLQGGI QLAQSRYWQT GDMYQGLGWE MLDWPVNPDS IINGS DNKIA LAARPVKAIT PPTPAVRASW VHKTGATGGF GSYVAFIPEK ELGIVMLANK NYPNPARVDA AWQILNALQ</pre>
予測される分子量	42 kDa including tags
領域	20 to 377
タグ	His tag N-Terminus
製品の詳細	Recombinant <i>E. coli</i> Beta-lactamase protein

特性

Our **Abpromise guarantee** covers the use of **ab104926** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

	Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid


前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 10% Glycerol (glycerin, glycerine), 0.316% Tris HCl
------------	---

関連情報

関連性	This protein is a serine beta-lactamase with a substrate specificity for cephalosporins.
-----	--

画像



15% SDS-PAGE analysis of 3ug ab104926

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors