

Recombinant Hantavirus Puumala, strain Vranica, nucleocapsid protein ab74555

[1 References](#) [画像数 1](#)

製品の詳細

製品名	Recombinant Hantavirus Puumala, strain Vranica, nucleocapsid protein
精製度	> 95 % SDS-PAGE. Purified by chromatography under denaturing conditions.
発現系	Yeast
アクセッション番号	A2ID68
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Hantavirus

特性

Our **Abpromise guarantee** covers the use of **ab74555** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA Western blot SDS-PAGE
製品の状態	Lyophilized

前処理および保存


保存方法および安定性	Shipped at 4°C. Store at +4°C. Constituent: PBS
再構成	Reconstitute with PBS or 8M urea to give 0.5mg/ml stock solution. After reconstitution in PBS aliquot and store at -20°C. Avoid freeze / thaw cycles. Proteins reconstituted in 8M urea could be stored at 4°C.

関連情報

関連性

Hantaviruses are rodent-borne agents that cause severe human diseases. The viruses constitute a group of antigenetically and genetically related envelope viruses with a tripartite, negative-sense RNA genome. In northern, western and central Europe the Puumala virus (PUUV) is a predominant serotype. It is the etiological agent of nephropathia epidemica, a mild form of hemorrhagic fever with renal syndrome.

画像



SDS-PAGE showing ab74555 at approximately 55kDa. (2 µg/lane).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors