

Recombinant HIV1 tat protein ab109929

製品の詳細

製品名	Recombinant HIV1 tat protein
精製度	> 95 % SDS-PAGE.
発現系	Escherichia coli
アクセッション番号	<u>Q6UFJ9</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
予測される分子量	11 kDa
領域	1 to 101

特性

Our **Abpromise guarantee** covers the use of **ab109929** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -80°C.
	Constituent: 0.1% Glycerol
再構成	The protein should be reconstituted in apyrogenic sterile water or PBS buffer. The reconstituted solution has to be used immediately. (Avoid repeated freeze/thaw cycles).

関連情報

関連性	HIV1 tat (Transactivator of Transcription) protein is a pleiotropic factor that induces a broad range of biological effects in numerous cell types. At the HIV promoter, tat is a powerful transactivator of gene expression which acts by inducing chromatin remodeling and by recruiting elongation-
-----	--

competent transcriptional complexes on to the viral LTR.

細胞内局在

Nuclear. Upon stimulation with EDN1, it is exported from the nucleus to the perinuclear region.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors