

Recombinant HSP65 protein ab92416

1 References

製品の詳細

製品名	Recombinant HSP65 protein
生理活性	This product does not contain GroEL.
精製度	> 90 % SDS-PAGE. ab92416 is purified by Multi-step chromatography
エンドトキシン・レベル	< 50.000 Eu/mg
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Mycobacterium bovis

特性

Our **Abpromise guarantee** covers the use of **ab92416** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot Functional Studies SDS-PAGE
製品の状態	Liquid
備考	This product does not contain GroEL.

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Constituents: 0.01999% Monobasic dihydrogen potassium phosphate, 0.11445% Dibasic monohydrogen sodium phosphate, 0.02025% Potassium chloride, 0.8004% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

関連情報

関連性

HSP65 is a member of the HSP60 family of heat shock proteins isolated from Mycobacterium bovis BCG. HSP65 from M. bovis is identical to that of M. tuberculosis, and similar to that of M. leprae, two highly pathogenic strains of mycobacterium. HSP65 is the immunodominant antigen during mycobacterial infection and vaccination, and has been linked to the development of autoimmune adjuvant arthritis in rats.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors