

Recombinant HSV1 gG Envelope Protein ab43048

1 References

製品の詳細

製品名	Recombinant HSV1 gG Envelope Protein
生理活性	Immunoreactive with sera of HSV-infected individuals, excellent antigen for detection of HSV with minimal specificity problems.
精製度	> 95 % SDS-PAGE. Sephrose derived purification.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
配列の追加情報	The protein contains the HSV-1 gG immunodominant regions.
製品の詳細	Recombinant HSV1 gG Envelope Protein

特性

Our **Abpromise guarantee** covers the use of **ab43048** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA Western blot
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. Constituents: 0.395% Tris HCl, 0.0292% EDTA, 50% Glycerol
------------	--

関連情報

関連性	Herpes simplex virus type 1 (HSV1) is usually associated with infections of the lips, mouth, and face. It is the most common herpes simplex virus and is usually acquired in childhood. HSV-1 often causes lesions inside the mouth such as cold sores (fever blisters) and is transmitted by
-----	---

contact with infected saliva. Glycoprotein G is suggested to contribute to viral entry through apical surfaces of polarized cells.

細胞内局在

Viral envelope

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors