

Recombinant HTLV I p24 protein ab43036

製品の詳細

製品名	Recombinant HTLV I p24 protein
生理活性	Immunoreactive with sera of HTLV I infected individuals, excellent antigen for early detection of HIV seroconvertors with minimal specificity problems.
精製度	> 95 % SDS-PAGE. Organic extraction > S-300 > HPLC.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
領域	314 to 344
配列の追加情報	The protein contains the full length sequence of HTLV-I p24.

特性

Our **Abpromise guarantee** covers the use of **ab43036** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA Western blot
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Store at -20°C. pH: 6.00 Constituents: 0.82% Sodium phosphate, 0.0154% DTT, 0.0292% EDTA
------------	--

関連情報

関連性	Human T-lymphotropic virus (HTLV) belongs to the family Retroviridae, genus Deltaretrovirus. It is a human, single-stranded RNA retrovirus that causes T-cell leukemia and T-cell lymphoma in adults and may also be involved in certain demyelinating diseases. Adult T-lymphotropic virus
-----	---

(ATLV) is a strain of this disease that effect primarily adults. Capsid protein p24 probably forms the conical core of the virus that encapsulates the genomic RNA-nucleocapsid complex.

細胞内局在

Core protein

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors