

Product datasheet

Recombinant Human 14-3-3 gamma/YWHAG protein ab53869

画像数 1

製品の詳細

製品名	Recombinant Human 14-3-3 gamma/YWHAG protein
精製度	> 95 % SDS-PAGE. Recombinant human 14-3-3 was expressed in E.coli and purified by using conventional chromatography techniques.
エンドトキシン・レベル	< 1.000 Eu/μg
発現系	Escherichia coli
アクセッション番号	<u>P61981</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MVDREQLVQK ARLAEQAERY DDMAAAMKNV TELNEPLSNE ERNLLSVAYK NVVGARRSSW RVISSIEQKT SADGNEKKIE MVRAYREKIE KELEAVCQDV LSLLDNYLIK NCSETQYESK VFYLMKMGDY YRYLAEVATG EKRAVVESS EKAYSEAHEI SKEHMQPHTP IRLGLALNYS VFYIEIQNAP EQACHLAKTA FDDAIAELDT LNEDSYKDST LIMQLLRDNL TLWTSDQDD DGEGNN
予測される分子量	28 kDa
領域	1 to 247

特性

Our **Abpromise guarantee** covers the use of **ab53869** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	This product was previously labelled as 14-3-3 gamma.

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 7.40

Constituents: PBS, 10% Glycerol (glycerin, glycerine)

関連情報

機能

Adapter protein implicated in the regulation of a large spectrum of both general and specialized signaling pathways. Binds to a large number of partners, usually by recognition of a phosphoserine or phosphothreonine motif. Binding generally results in the modulation of the activity of the binding partner.

組織特異性

Highly expressed in brain, skeletal muscle, and heart.

配列類似性

Belongs to the 14-3-3 family.


翻訳後修飾

Phosphorylated by various PKC isozymes.

細胞内局在

Cytoplasm.

画像



3µg by SDS-PAGE under reducing conditions and visualized by coomassie blue stain.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you

- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors