

Recombinant Human ADRM1/ARM-1 protein ab124570

画像数 1

製品の詳細

製品名	Recombinant Human ADRM1/ARM-1 protein
精製度	> 90 % SDS-PAGE. ab124570 was purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	Q16186
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHSSGLVPRGSHMGSHMTTSGALFPSLVPGSR GASNKYLVEF RAGKMSLKGTTPDKRKGVLVYIQQTDDSLIHFCWKDRSGN VEDDLIIF PDDCEFKRVPQCPSGRVYVLKFKAGSKRLLFFWMQEPKTDQDE EHCRCVNE YLNPPMPGALGASGSSGHELSALGGEGGLQSLLGNMHSQ MQLIGPAG LGGLGGLGALTGPGLASLLGSSGPPGSSSSSSRSQSAAVT SSTTSSTR ATPAPSAPAAASATSPSPAPSSGNGASTAASPTQPIQLSDLQ SILATMNV PAGPAGGQQVDLASVLTPEIMAPILANADVQERLLPYLP SLPQTAE IQNTLTSPQFQALGMFSAALASGQLGPLMCQFGLPAEAVEA ANKGDVEA FAKAMQNNAKPEQKEGDTKDKKDEEEDMSLD
予測される分子量	45 kDa including tags
領域	1 to 407
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab124570** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid
備考	This product was previously labelled as ADRM1


前処理および保存

保存方法および安定性	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. pH: 8.00 Constituents: 0.02% DTT, 0.32% Tris HCl, 20% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
------------	--

関連情報

機能	Functions as a proteasomal ubiquitin receptor. Recruits the deubiquitinating enzyme UCHL5 at the 26S proteasome and promotes its activity.
配列類似性	Belongs to the ADRM1 family. Contains 1 PH domain.
ドメイン	The PH domain mediates interactions with PSMD1 and ubiquitin. Preferential binding to the proximal subunit of K48-linked diubiquitin allows UCHL5 access to the distal subunit.
細胞内局在	Cytoplasm. Nucleus.

画像



15% SDS-PAGE analysis of 3 µg ab124570.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors