

Recombinant Human ADX protein ab87670

画像数 1

製品の詳細

製品名	Recombinant Human ADX protein
精製度	> 90 % SDS-PAGE. ab87670 is purified using conventional chromatography techniques. Endotoxin: < 1.0 EU per 1 µg of protein (determined by LAL method).
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MASMTGGGQM GRGSMSSSED KITVHFINRD GETLTTKGKV GDSLLDVVVE NNLDIDGFGA CEGTLACSTC HLIFEDHIYE KLDAITDEEN DMLDLAYGLT DRSRLGCQIC LTKSMDNMTV RVPETVADAR QSIDVGKTS

特性

Our **Abpromise guarantee** covers the use of **ab87670** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	This product was previously labelled as Adrenodoxin


前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)
------------	--

関連情報

機能	Participates in the synthesis of thyroid hormones. Transfers electrons from adrenodoxin reductase to the cholesterol side chain cleavage cytochrome P450.
配列類似性	Belongs to the adrenodoxin/putidaredoxin family. Contains 1 2Fe-2S ferredoxin-type domain.
細胞内局在	Mitochondrion matrix.

画像



15% SDS-PAGE showing ab87670 at approximately 15kDa (3μg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors