

Recombinant Human AF6 protein ab91389

製品の詳細

製品名	Recombinant Human AF6 protein
精製度	> 95 % SDS-PAGE.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
領域	1 to 141

特性

Our **Abpromise guarantee** covers the use of **ab91389** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. pH: 7.50 Constituents: 0.595% HEPES, 0.0292% EDTA, 50% Glycerol, 0.29% Sodium chloride
------------	---

関連情報

機能	Belongs to an adhesion system, probably together with the E-cadherin-catenin system, which plays a role in the organization of homotypic, interneuronal and heterotypic cell-cell adherens junctions (AJs). Nectin- and actin-filament-binding protein that connects nectin to the actin cytoskeleton.
関連疾患	Note=A chromosomal aberration involving MLLT4 is associated with acute leukemias. Translocation t(6;11)(q27;q23) with MLL/HRX. The result is a rogue activator protein.

配列類似性	Contains 1 dilute domain. Contains 1 FHA domain. Contains 1 PDZ (DHR) domain. Contains 2 Ras-associating domains.
ドメイン	The PDZ/DHR domain interacts with the C-terminus of nectin and the Pro-rich N-terminus domain interacts with F-actin.
細胞内局在	Cell junction > adherens junction. Not found at cell-matrix AJs.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors