

# Recombinant Human Alpha SNAP protein ab100963

画像数 1

### 製品の詳細

製品名	Recombinant Human Alpha SNAP protein
精製度	> 95 % SDS-PAGE. ab100963 is purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<b><u>P54920</u></b>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<p>MGSSHHHHHHSSGLVPRGSHMDNSGKEAEAMALLAEAERKVK  NSQSFFSG  LFGSSKIEEACEIYARAANMFKMAKNWSAAGNAFCQAAQLH  LQLQSKHD  AATCFVDAGNAFKKADPQEAINCLMRAIEIYTDAGRFTIAAK  HHISIAEI  YETELVDIEKAIAHYEQSADYYKGEESSANKCLLKVAGYA  ALLEQYQK  AIDIYEQVGTNAMDSPLLKYSADYFFKAALCHFCIDMLNAK  LAVQKYEE  LFPAFSDSRECKLMKKLLEAHEEQNVDSYTESVKEYDSISRL  DQWLTTML LRIKKTIQGDEEDLR</p>
予測される分子量	35 kDa including tags
領域	1 to 295
タグ	His tag N-Terminus

### 特性

Our **Abpromise guarantee** covers the use of **ab100963** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Mass Spectrometry
質量分析	MALDI-TOF

## 製品の状態

Liquid

## 前処理および保存

### 保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

pH: 7.50

Constituents: 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)

## 関連情報

### 機能

Required for vesicular transport between the endoplasmic reticulum and the Golgi apparatus.


### 配列類似性

Belongs to the SNAP family.

### 細胞内局在

Membrane.

## 画像



15% SDS-PAGE showing ab100963 (3 µg) at approximately 35.3 kDa.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- 
- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors