

# Recombinant Human Annexin-6/ANXA6 protein ab92968

画像数 2

### 製品の詳細

---

製品名	Recombinant Human Annexin-6/ANXA6 protein
精製度	> 95 % SDS-PAGE. Purified by 6xHis tag / Immobilized-metal affinity chromatography (IMAC).
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
領域	1 to 673

### 特性

---

Our **Abpromise guarantee** covers the use of **ab92968** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot ELISA SDS-PAGE Mass Spectrometry
製品の状態	Liquid
備考	This product was previously labelled as Annexin VI

### 前処理および保存


---

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Preservative: 0.68% Imidazole Constituents: 0.0078% Beta mercaptoethanol, 0.316% Tris HCl, 2.9% Sodium chloride
------------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## 関連情報

機能	May associate with CD21. May regulate the release of Ca(2+) from intracellular stores.
配列類似性	Belongs to the annexin family. Contains 8 annexin repeats.
ドメイン	A pair of annexin repeats may form one binding site for calcium and phospholipid.
翻訳後修飾	Phosphorylated in response to growth factor stimulation.
細胞内局在	Cytoplasm. Melanosome. Identified by mass spectrometry in melanosome fractions from stage I to stage IV.

## 画像




10% SDS-PAGE showing ab92968

Lane 1 MWM

Lane 2 Stained with Coomassie Blue

Lane 3 Immunoblotting with anti-6xHis serum and peptide fingerprinting by MALDI-TOF-TOF mass spectrometry

SDS-PAGE - Recombinant Human Annexin-6/ANXA6 protein (ab92968)



Western blot - Recombinant Human Annexin-6/ANXA6 protein (ab92968)

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors