

Recombinant Human Annexin-8/ANXA8 protein ab99215

画像数 1

製品の詳細

製品名	Recombinant Human Annexin-8/ANXA8 protein
精製度	> 90 % SDS-PAGE. ab99215 was purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<u>P13928</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<p>MGSSHHHHHSSGLVPRGSHMAWWKSWIEQEGVTVKSSSH FNPDPDAETL YKAMKGIGTNEQAIIDVLTKRSNTQRQQIAKSFKAQFGKDLT ETLKSELS GKFERLIVALMYPPRYEAKELHDAMKGLGTKEGVIIIEILAS RTKNQLRE IMKAYEEDYGSSLEEDIQADTSGYLERILVCLLQGSRDDVSS FVDPGLAL QDAQDLYAAGEKIRGTDEMKFITILCTRSATHLLRVFEEYEK IANKSIED SIKSETHGSLEEAMLTVVKCTQNLHSYFAERLYYAMKGAGTR DGTLLIRNI VSRSEIDLNLKCHFKKMYGKTLSSMIMEDTSGDYKNALLSL VGSDP</p>
予測される分子量	39 kDa including tags
領域	1 to 327
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab99215** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション Mass Spectrometry

	SDS-PAGE
質量分析	MALDI-TOF
製品の状態	Liquid
備考	This product was previously labelled as Annexin VIII


前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.0154% DTT, 0.316% Tris HCl, 0.0292% EDTA, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride
------------	--

関連情報

機能	This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade.
配列類似性	Belongs to the annexin family. Contains 4 annexin repeats.
ドメイン	A pair of annexin repeats may form one binding site for calcium and phospholipid.

画像



15% SDS-PAGE analysis of 3µg ab99215.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors