

Recombinant human Axl protein ab64298

画像数 5

製品の詳細

製品名	Recombinant human Axl protein
生理活性	Specific activity: 104 nmol/min/mg.
精製度	> 90 % SDS-PAGE. Affinity purified.
発現系	Baculovirus infected Sf9 cells
アクセッション番号	<u>P30530</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
領域	473 to 894

特性

Our **Abpromise guarantee** covers the use of **ab64298** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Functional Studies
製品の状態	Liquid
備考	<u>ab204853</u> (IRS1 peptide) can be utilized as a substrate for assessing kinase activity


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7.00 Preservative: 1.02% Imidazole Constituents: 0.00174% PMSF, 0.82% Sodium phosphate, 0.00308% DTT, 25% Glycerol (glycerin, glycerine), 1.74% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	--

関連情報


機能	May function as a signal transducer between specific cell types of mesodermal origin. In case of filovirus infection, seems to function as a cell entry factor.
組織特異性	Highly expressed in metastatic colon tumors. Expressed in primary colon tumors. Weakly expressed in normal colon tissue.
配列類似性	Belongs to the protein kinase superfamily. Tyr protein kinase family. AXL/UFO subfamily. Contains 2 fibronectin type-III domains. Contains 2 Ig-like C2-type (immunoglobulin-like) domains. Contains 1 protein kinase domain.
細胞内局在	Membrane.

画像




The specific activity of Axl (ab64298) was determined to be 88 nmol/min/mg as per activity assay protocol

Functional Studies - Recombinant human Axl protein (ab64298)





SDS PAGE analysis of ab64298

SDS-PAGE - Recombinant human Axl protein (ab64298)



SDS PAGE analysis of ab64298



SDS Page analysis of ab64298.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors