

Recombinant Human BAALC protein ab95394

画像数 1

製品の詳細

製品名	Recombinant Human BAALC protein
精製度	> 95 % SDS-PAGE. ab95394 is purified using conventional chromatography techniques.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHH SSGLVPRGSH MCGGSRADA IEPRYYESWT RETESTWLTY TDSDAPPSAA APDSGPEAGG LHSGMLEDGL PSNGVPRSTA PGGIPNPEKK TNCETQCPNP QSLSSGPLTQ KQNGLQTTEA KRDAKRMPAK EVTINVTDSI QQMDRSRRIT KNCVN
領域	1 to 145

特性

Our **Abpromise guarantee** covers the use of **ab95394** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Mass Spectrometry
	SDS-PAGE

製品の状態	Liquid
-------	--------


前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. pH: 8.00 Constituents: 0.0154% DTT, 0.316% Tris HCl, 0.0584% EDTA, 20% Glycerol (glycerin, glycerine)
------------	---

関連情報

機能	May play a synaptic role at the postsynaptic lipid rafts by interacting with CAMK2A.
組織特異性	Expressed by neural and hematopoietic cells. Expression is found both in normal early progenitor cells and in the most immature type of blasts in acute leukemia but not in mature hematopoietic cells. Isoform 2 and isoform 6 are expressed in the brain. Isoform 2 shows a low expression in neuroectoderm-derived tissues such as adrenal gland and no expression in bone marrow, peripheral blood lymphocytes or lymph nodes, or in tumors or cancer cell lines of nonneural tissue origin.
翻訳後修飾	Palmitoylation and myristoylation target the protein to the lipid rafts.
細胞内局在	Cytoplasm. Membrane. Postsynaptic lipid rafts.

画像



15% SDS-PAGE showing ab95394 at approximately 17.7kDa (3μg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors