

Recombinant Human Bak protein ab114337

画像数 1

製品の詳細

製品名	Recombinant Human Bak protein
発現系	Wheat germ
アクセッション番号	<u>Q16611</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MASGQGPGRPRQECGEPALPSASEEQVAQDTEEVFRSYVFYR HQQEAE GVAAPADPEMVTLPQLPSSTMGQVGRQLAIIGDDINRRYDSE FQTMLQHL QPTAENAYEYFTKIATSLFESGINWGRVVALLGFGYRLALHV YQHGLTGF LGQVTRFVDFMLHHCIARWIAQRGGWVAALNLGNPILNVL VVLGVLL GQFVRRFFKS
予測される分子量	49 kDa including tags
領域	1 to 211
タグ	GST tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab114337** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE
	ELISA
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
------------	---


pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

関連情報

機能	In the presence of an appropriate stimulus, accelerates programmed cell death by binding to, and antagonizing the anti-apoptotic action of BCL2 or its adenovirus homolog E1B 19k protein. Low micromolar levels of zinc ions inhibit the promotion of apoptosis.
組織特異性	Expressed in a wide variety of tissues, with highest levels in the heart and skeletal muscle.
配列類似性	Belongs to the Bcl-2 family.
ドメイン	Intact BH3 motif is required by BIK, BID, BAK, BAD and BAX for their pro-apoptotic activity and for their interaction with anti-apoptotic members of the Bcl-2 family.
細胞内局在	Mitochondrion membrane.

画像



ab114337 analysed on a 12.5% SDS-PAGE gel stained with Coomassie Blue

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

-
- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors