

Recombinant Human BIT1 protein ab101645

1 References [画像数 1](#)

製品の詳細

製品名	Recombinant Human BIT1 protein
精製度	> 95 % SDS-PAGE. ab101645 was purified using anion-exchange chromatography (DEAE sepharose resin) and gel-filtration chromatography (Sephacryl S-200) with 20mM Tris pH 7.5, 2mM EDTA.
発現系	Escherichia coli
アクセッション番号	<u>Q9Y3E5</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHSSGLVPRGSHMEYKMILVVRNDLKMGGKGV AAQCSHAAVS AYKQIQRRNPEMLKQWEYCGQPKVVVKAPDEETLIALAHAK MLGLTVSL IQDAGRTQIAPGSQTVLGIGPGPADLIDKVTGHLKLY
予測される分子量	15 kDa including tags
領域	64 to 179
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab101645** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid
備考	This product was previously labelled as PTRH2

前処理および保存

保存方法および安定性

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 0.0154% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)

関連情報

機能

The natural substrate for this enzyme may be peptidyl-tRNAs which drop off the ribosome during protein synthesis.

Promotes caspase-independent apoptosis by regulating the function of two transcriptional regulators, AES and TLE1.


配列類似性

Belongs to the PTH2 family.

細胞内局在

Mitochondrion.

画像



ab101645 at 3 µg analysed by 15% SDS PAGE

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors