

Product datasheet

Recombinant Human BNIP3L/NIX protein ab114290

画像数 1

製品の詳細

製品名	Recombinant Human BNIP3L/NIX protein
発現系	Wheat germ
アクセッション番号	O60238
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	SSNGNDNGNGKNGGLEHVPSSSSIHNGDMEKILLDAQHESGQ SSSRGSSH CDSPSPQEDGQIMFDVEMHTSRDHSSQSEEEVVEGEKE
予測される分子量	35 kDa including tags
領域	43 to 130

特性

Our **Abpromise guarantee** covers the use of **ab114290** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot ELISA SDS-PAGE
製品の状態	Liquid
備考	This product was previously labelled as BNIP3L.


前処理および保存

保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
------------	---

関連情報

機能	Induces apoptosis. Interacts with viral and cellular anti-apoptosis proteins. Can overcome the suppressors BCL-2 and BCL-XL, although high levels of BCL-XL expression will inhibit apoptosis. May function as a tumor suppressor. Inhibits apoptosis induced by BNIP3.
配列類似性	Belongs to the NIP3 family.
翻訳後修飾	Undergoes progressive proteolysis to an 11 kDa C-terminal fragment, which is blocked by the proteasome inhibitor lactacystin.
細胞内局在	Nucleus envelope. Endoplasmic reticulum. Mitochondrion. Membrane.

画像



ab114290 analysed on a 12.5% SDS-PAGE gel stained with Coomassie Blue.

SDS-PAGE - Recombinant Human BNIP3L/NIX protein (ab114290)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors