

Recombinant Human Calcipressin 1/RCAN1 protein ab63837

画像数 1

製品の詳細

製品名	Recombinant Human Calcipressin 1/RCAN1 protein
精製度	> 90 % SDS-PAGE. ab63837 was expressed in E.coli and purified by using conventional chromatography techniques.
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	MKLYFAQTLH IGSSHLAPPN PDKQFLISPP ASPPVGWKQV EDATPVINYD LLYAISKLG GEKYELHAAT DTTSPVVVHV CESDQEKEE EEMERMRRPK PKIIQTRRPE YTPIHLS
領域	1 to 117

特性

Our **Abpromise guarantee** covers the use of **ab63837** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	This product was previously labelled as Calcipressin 1


前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 7.5 Constituents: 0.0154% DTT, 0.242% Tris, 0.58% Sodium chloride
------------	---

関連情報

機能	Inhibits calcineurin-dependent transcriptional responses by binding to the catalytic domain of calcineurin A. Could play a role during central nervous system development.
組織特異性	Highly expressed heart, brain and skeletal muscle. Also expressed in all other tissues.
配列類似性	Belongs to the RCAN family.

画像



SDS Page analysis of ab63837 (3 ug on 15% gel).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors