

Product datasheet

Recombinant Human CARHSP1 protein ab101171

画像数 1

製品の詳細

製品名	Recombinant Human CARHSP1 protein
精製度	> 90 % SDS-PAGE. ab101171 is purified using conventional chromatography techniques.
発現系	Escherichia coli
アクセッション番号	<u>Q9Y2V2</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHH SSGLVPRGSH MSSEPPPPPQ PPTHQASVGL LDTPRSRERS PSPLRGNVVP SPLPTRRTRT FSATVRASQG PVYKGVCKCF CRSKGHGFIT PADGGPDIFL HISDVEGEYV PVEGDEVYK MCSIPPKNEK LQAVEVVITH LAPGTKHETW SGHVISS
予測される分子量	18 kDa
領域	1 to 147
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab101171** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid

前処理および保存

保存方法および安定性 Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -

80°C. Avoid freeze / thaw cycle.


pH: 8.00

Constituents: 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.58% Sodium chloride

関連情報

機能	Binds mRNA and regulates the stability of target mRNA. Binds single-stranded DNA (in vitro).
配列類似性	Contains 1 CSD (cold-shock) domain.
翻訳後修飾	Dephosphorylated by calcineurin in a Ca(2+) dependent manner.
細胞内局在	Cytoplasm. Cytoplasm > P-body. Cytoplasmic granule. Detected at cytoplasmic stress granules and P-bodies. Detected at exosome granules where mRNA is degraded.

画像



15% SDS-PAGE analysis of ab101171 (3ug)

SDS-PAGE - Recombinant Human CARHSP1
protein (ab101171)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors