

Recombinant Human CCL1 protein ab104828

画像数 1

製品の詳細

製品名	Recombinant Human CCL1 protein
精製度	> 85 % SDS-PAGE. ab104828 is purified using conventional chromatography techniques, after refolding of the isolated inclusion bodies in a renaturation buffer.
発現系	Escherichia coli
アクセッション番号	<u>P22362</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHHSSGLVPRGSHMKSMQVPFSRCCFSFAEQEI PLRAILCYRN TSSICSNEGLIFKLKRGKEACALDTVGVQRRKMLRHCP SKRK
予測される分子量	11 kDa including tags
領域	24 to 96
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab104828** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	Mass Spectrometry
質量分析	MALDI-TOF
製品の状態	Liquid
備考	This product was previously labelled as I-309

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


pH: 7.50

Constituents: 0.0308% DTT, 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine), 0.29% Sodium chloride

関連情報

機能	Cytokine that is chemotactic for monocytes but not for neutrophils. Binds to CCR8.
配列類似性	Belongs to the intercrine beta (chemokine CC) family.
細胞内局在	Secreted.

画像



15% SDS-PAGE showing ab104828 at approximately 10.8kDa (3μg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors