

Recombinant human CCL27 protein ab50060

製品の詳細

製品名	Recombinant human CCL27 protein
生理活性	Determined by its ability to chemoattract CXCR3 transfected cells using a concentration of 10.0 - 100.0 ng/mL.
精製度	> 98 % SDS-PAGE. Greater than 98% by SDS-PAGE gel and HPLC analyses.
エンドトキシン・レベル	< 0.100 Eu/μg
発現系	Escherichia coli
アクセッション番号	<u>Q9Y4X3</u>
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	FLLPPSTACC TQLYRKPLSD KLLRKVIQVE LQEADGDCHL QAFVLHLAQR SICIHPQNPS LSQWFEHQER KLHGTLPKLN FGMLRKMG
予測される分子量	10 kDa
領域	25 to 112
配列の追加情報	Full length mature chain without signal peptide.

特性

Our **Abpromise guarantee** covers the use of **ab50060** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Functional Studies SDS-PAGE
製品の状態	Lyophilized

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	---

再構成

Reconstitute to 1 mg/ml with distilled water.

関連情報

機能	Chemotactic factor that attracts skin-associated memory T-lymphocytes. May play a role in mediating homing of lymphocytes to cutaneous sites. Binds to CCR10.
組織特異性	Testis, thymus, placenta, ovary and skin.
配列類似性	Belongs to the intercrine beta (chemokine CC) family.
細胞内局在	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors