

Recombinant Human CD1d protein ab114936

画像数 1

製品の詳細

製品名	Recombinant Human CD1d protein
発現系	Wheat germ
アクセッション番号	P15813
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGCLLFLLWALLQAWGSAEVPQRLFPLRCLQISSFANSSWT RTDGLAWL GELQTHSWSNDSDTVRSCLKPWSQGTFFSDQQWETLQHIFRVYR SSFTRDVK EFAKMLRLSYPLELQVSAGCEVHPGNASNNFFHVAFAQKDIL SFQGSWE PTQEAPLWVNLAIQVLNQDKWTRETVQWLLNGTCPQFVSGLL ESGKSELK KQVKPKAWLSRGPSPGPGRLLLCHVSGFYKPVVVKWMRGE QEQQGTQP GDILPNADETWYLRATLDVVAGEAAGLSCRVKHSSLEGQDIV LYWGGSYT SMGLIALAVLACLFLLLIVGFTSRFKRQTSYQGVL
予測される分子量	63 kDa including tags
領域	1 to 335
タグ	GST tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab114936** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Western blot
	SDS-PAGE
	ELISA
製品の状態	Liquid

備考 Protein concentration is lot specific and is between 0.01 - 0.12 ug/ul

前処理および保存

保存方法および安定性 Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
pH: 8.00
Constituents: 0.3% Glutathione, 0.79% Tris HCl

関連情報


機能 Antigen-presenting protein that binds self and non-self glycolipids and presents them to T-cell receptors on natural killer T-cells.

組織特異性 Expressed on cortical thymocytes, on certain T-cell leukemias, and in various other tissues.

配列類似性 Contains 1 Ig-like (immunoglobulin-like) domain.

細胞内局在 Cell membrane. Endosome membrane. Lysosome membrane. Subject to intracellular trafficking between the cell membrane, endosomes and lysosomes.

画像



ab114936 on a 12.5% SDS-PAGE Stained with Coomassie Blue.

SDS-PAGE - Recombinant Human CD1d protein
(ab114936)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors