

Recombinant Human CD70 protein ab119815

製品の詳細

製品名	Recombinant Human CD70 protein
精製度	> 95 % SDS-PAGE. Purity: Greater than 95% by SDS-PAGE gel and HPLC analyses
エンドトキシン・レベル	< 0.100 Eu/μg
発現系	CHO cells
アクセッション番号	<u>P32970</u>
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	HHHHHHHPS PGGSGGQRFA QAQQQLPLES LGWDVAELQL NHTGPQQDPR LYWQGGPALG RSFLHGPELD KGQLRIHRDG IYMVHIQVTL AICSSTTASR HHPTTLAVGI CSPASRSISL LRLSFHQGCT IASQRLTPLA RGDTLCTNLT GTLLPSRNTD ETFFGVQWVR P
領域	39 to 193
タグ	His tag N-Terminus

特性

Our **Abpromise guarantee** covers the use of **ab119815** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
	HPLC
	Functional Studies

製品の状態	Lyophilized
-------	-------------

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
------------	---

Constituent: 0.6% Sodium citrate

再構成

For lot specific reconstitution information please contact our Scientific Support Team.

関連情報

機能	Cytokine that binds to CD27. Plays a role in T-cell activation. Induces the proliferation of costimulated T-cells and enhances the generation of cytolytic T-cells.
配列類似性	Belongs to the tumor necrosis factor family.
細胞内局在	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors