

Recombinant Human CDX2 protein ab114247

画像数 2

製品の詳細

製品名	Recombinant Human CDX2 protein
発現系	Wheat germ
アクセッション番号	Q99626
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	<p>MYVSYLLDKDVSMYPSSVRHSGGLNLAPQNFVSPQYPDYGG YHVAAAAA AAANLDSAQSPGPSWPAAYGAPLREDWNGYAPGGAAAAANAV AHGPNNGGS PAAAMGYSSPADYHPHHHPHHHPHHPAAAPSCASGLLQTLNP GPPGPAAT AAAELSPGGQRRNLCEWMRKPAQQLGSQVKTRTKDKYRVV YTDHQRLE LEKEFHYSRYITIRRKAEALATLGLSERQVKIWFQNRRAKER KINKKKLQ QQQQQQPPPPPPPPPPQPPQPGPLRSVPEPLSPVSSLQASV SGSVPGVL GPTGGVLNPTVTQ</p>
予測される分子量	61 kDa including tags
領域	1 to 313

特性

Our **Abpromise guarantee** covers the use of **ab114247** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	ELISA
	SDS-PAGE
	Western blot
製品の状態	Liquid

備考

前処理および保存

保存方法および安定性

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

関連情報

機能

Involved in the transcriptional regulation of multiple genes expressed in the intestinal epithelium. Important in broad range of functions from early differentiation to maintenance of the intestinal epithelial lining of both the small and large intestine.

配列類似性

Belongs to the Caudal homeobox family.

Contains 1 homeobox DNA-binding domain.


翻訳後修飾

Phosphorylation of Ser-60 mediates the transactivation capacity.

細胞内局在


Nucleus.

画像



SDS-PAGE - Recombinant Human CDX2 protein
(ab114247)

12.5% SDS-PAGE showing ab114247 at approximately 60.54kDa stained with Coomassie Blue.



Anti-CDX2 antibody (**ab88129**) at 1 µg/ml + Recombinant Human CDX2 protein (ab114247) at 0.01 µg

Secondary

Goat Anti-Rabbit IgG H&L (HRP) preadsorbed (**ab97080**) at 1/5000 dilution

Developed using the ECL technique.

Performed under reducing conditions.

Exposure time: 10 seconds

ab88129 recognizes the full length tagged recombinant CDX2 protein (ab114247) which has an expected molecular weight of 61 kDa.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors