

Recombinant human Chk1 protein ab69918

画像数 5

製品の詳細

製品名	Recombinant human Chk1 protein
生理活性	Specific activity: 144 nmol /min/mg as determined by Kinase assay.
精製度	> 80 % Densitometry. Affinity purified.
発現系	Baculovirus infected Sf9 cells
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human

特性

Our **Abpromise guarantee** covers the use of **ab69918** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE Functional Studies
製品の状態	Liquid
備考	ab204854 (Cdc25C peptide) can be utilized as a substrate for assessing kinase activity

前処理および保存


保存方法および安定性	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 7 Constituents: 0.0038% EGTA, 0.00174% PMSF, 0.00385% DTT, 0.79% Tris HCl, 0.00292% EDTA, 25% Glycerol (glycerin, glycerine), 0.87% Sodium chloride This product is an active protein and may elicit a biological response in vivo, handle with caution.
------------	---

関連情報



The specific activity of Chk1 (ab69918) was determined to be 165 nmol/min/mg as per activity assay protocol


Functional Studies - Recombinant human Chk1 protein (ab69918)




SDS PAGE analysis of ab69918

SDS-PAGE - Recombinant human Chk1 protein (ab69918)


SDS PAGE analysis of ab69918



SDS-PAGE - Recombinant human Chk1 protein
(ab69918)



Functional Studies - Recombinant human Chk1
protein (ab69918)



SDS-PAGE - Recombinant human Chk1 protein
(ab69918)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors