

Recombinant Human Chromogranin A protein ab85339

画像数 1

製品の詳細

製品名	Recombinant Human Chromogranin A protein
精製度	> 95 % SDS-PAGE. ab85339 is purified using conventional chromatography techniques. Endotoxin Level: < 1.0 EU per 1 µg of protein (determined by LAL method).
発現系	Escherichia coli
タンパク質長	Protein fragment
Animal free	No
由来	Recombinant
生物種	Human
配列	MLPVNSPMNK GDTEVMKCIV EVISDTLSKP SPMPVSQECF ETLRGDERIL SILRHQNLLK ELQDLALQGA KERAHQQKKH SGFEDELSEV LENQSSQAEL KEAVEEPSSK DVME
領域	19 to 131
配列の追加情報	This fragment is from the N terminus of Chromogranin A.

特性

Our **Abpromise guarantee** covers the use of **ab85339** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション SDS-PAGE

製品の状態 Liquid

備考

Vasostatin 2 is the N-terminal fragment (19-131 aa) derived from the cleavage of chromogranin A (CgA) which is a member of the chromogranin/secretogranin(granins) family of neuroendocrine secretory proteins. Vasostatin 2 has been shown to exert several biological activities on several tissues and organs and exerts a large spectrum of homeostatic actions, including antifungal and antimicrobial effect, modulation of cell adhesion, and inhibition of parathyroid hormone secretion.

前処理および保存

保存方法および安定性

Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.


pH: 8.00

Constituents: 0.316% Tris HCl, 10% Glycerol (glycerin, glycerine)

関連情報

機能	Pancreastatin strongly inhibits glucose induced insulin release from the pancreas.
配列類似性	Belongs to the chromogranin/secretogranin protein family.
翻訳後修飾	Sulfated on tyrosine residues and/or contains sulfated glycans. O-glycosylated with core 1 or possibly core 8 glycans.
細胞内局在	Secreted. Neuroendocrine and endocrine secretory granules.

画像



15% SDS-PAGE showing ab85339 at approximately 12.8kDa (3μg).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors