

Recombinant Human CIB1/KIP protein ab78759

画像数 1

製品の詳細

製品名	Recombinant Human CIB1/KIP protein
精製度	> 95 % SDS-PAGE. ab78759 is purified by using anion-exchange chromatography (DEAE sepharose resin) and gel-filtration chromatography (Sephacryl S-200) with 20mM Tris pH 7.5, 2mM EDTA.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHH SSGLVPRGSH MGGSGSRLSK ELLAEYQDLT FLTKQEILLA HRRFCELLPQ EQRTVESSLR AQVPFEQILS LPELKANPFK ERICRVFSTS PAKDSLSFED FLDLLSVFSD TATPDIKSHY AFRIFD FDDD GTLNREDLSR LVNCLTGEGE DTRLSASEMK QLIDNILEES DIDRDGTINL SEFQHVISRS PDFASSFKIV L

特性

Our **Abpromise guarantee** covers the use of **ab78759** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid
備考	This product was previously labelled as CIB1


前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.077% DTT, 0.242% Tris
------------	--

関連情報

機能	May convert the inactive conformation of integrin α -IIb/ β 3 to an active form through binding to the integrin cytoplasmic domain. Induces cell migration and spreading mediated through integrin (possibly via focal adhesion complexes). Functions as a negative regulator of stress activated MAP kinase (MAPK) signaling pathways. May play a role in regulation of apoptosis. Interacts with and up-regulates PTK2 activity. Down regulates inositol 1,4,5-trisphosphate receptor-dependent calcium signaling.
組織特異性	Ubiquitous.
配列類似性	Contains 2 EF-hand domains.
細胞内局在	Membrane. Cytoplasm. Nucleus. Cell projection > filopodium. Apical cell membrane. Localizes to the perinuclear region in the presence of NBR1. Colocalizes with TAS1R2 in apical regions of taste receptor cells.

画像



15% SDS-PAGE showing ab78759 at approximately 23kDa (3 μ g).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors