

Recombinant Human CPI-17 protein ab75206

画像数 1

製品の詳細

製品名	Recombinant Human CPI-17 protein
精製度	> 95 % SDS-PAGE. ab75206 is purified by using conventional chromatography techniques.
発現系	Escherichia coli
タンパク質長	Full length protein
Animal free	No
由来	Recombinant
生物種	Human
配列	MGSSHHHHHH SSGLVPRGSH MAAQRLGKRV LSKLQSPSRA RGPGGSPGGL QKRHARVTVK YDRRELQRRRL DVEKWIDGRL EELYRGMEAD MPDEINIDEL LELESEEEERS RKIQGLLKSC GKPVEDFIQE LLAQLQGLHR QPGLRQPSPS HDGSLSPQLQD RARTAHP

特性

Our **Abpromise guarantee** covers the use of **ab75206** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.0154% DTT, 0.242% Tris, 0.0584% EDTA, 10% Glycerol (glycerin, glycerine)
------------	---

関連情報

機能	Inhibitor of PPP1CA. Has over 1000-fold higher inhibitory activity when phosphorylated, creating a
----	--

組織特異性

molecular switch for regulating the phosphorylation status of PPP1CA substrates and smooth muscle contraction.


配列類似性

Isoform 1 is detected in aorta and testis. Isoform 2 is detected in aorta.

細胞内局在

Belongs to the PP1 inhibitor family.

Cytoplasm.

画像

15% SDS-PAGE showing ab75206 at approximately 18kDa. (3μg)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors